

pH 4, pH 7 & pH 10

Buffer Capsules (100ml pH Buffer Standard)

The most comprehensive line of pH Buffer standards used for laboratory and analytical measurement with 101 buffer values.

- Hydrion Buffer salts are stable, and easy to prepare. Accurate to known pH values at ± 0.02 @ 25°C.

Preparation:

To prepare the Buffer Standard, empty the contents of the Hydrion pH Buffer Capsule into 100ml de-ionized or distilled water.

Storage And Shelf Life:

When stored at room temperature, in clean airtight bottles, prepared solutions will last for three weeks. When added to the buffer solution, our Hydrion Buffer Preservative will protect prepared solutions for up to three months.

Unprepared, packaged Hydrion Buffer Salts have a shelf life of five years.

- All HYDRION BUFFER SALTS ARE DIRECTLY TRACEABLE TO NIST STANDARDS AND CERTIFIED TO ± 0.02 pH UNITS AT 25°C



pH 4.00 Buffer, 10 Capsule of vial



pH 7.00 Buffer, 10 Capsule of vial



pH 10.00 Buffer, 10 Capsule of vial

Ordering

pH 4 : 4122F21.01 (10 capsule)

pH 7 : 4122F36.01 (10 capsule)

pH10 : 4122F51.01 (10 capsule)

Website : www.eastco-hk.com

EASTCO LIMITED

儀 高 儀 器 有 限 公 司



香港總公司
香港新界荃灣沙咀道26-38號
匯力工業中心2樓26室
Tel : (852)2456 1823
Fax: (852) 2469 1767
E-mail: info@eastco-hk.com

中國公司
儀高南儀器(深圳)有限公司
中國深圳市八卦四路泰克大廈
3樓V區 郵編: 518029
Tel : (86-755) 8241 3890
Fax: (86-755) 8241 3857
E-mail: shenzhen@eastco-hk.com

上海辦事處
中國上海市滬太路701號
18樓1810室
郵編: 200072
Tel : (86-21) 5653 6311
Fax: (86-21) 5653 6309
E-mail: shanghai@eastco-hk.com